## Upute za uporabu

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### PODACI ZA KONTAKT

Jiaxing Funing Medical Instrument Co.,Ltd.  
3<sup>rd</sup> Floor, Building C, Innovation town, Guigu Park,  
Luoxing Street, 314100 Jiashan County, Zhejiang Province, P.R. Kina

Prolinx GmbH  
Brehmstr. 56, 40239, Duesseldorf  
Tel.: 0049 2131 4051968-0  
Faks: 0049 2131 4051968-9  
Osoba za kontakt: G. Nianzhua Liu  
Dimdi kod: DE/0000045300  
E-pošta: dusai@eulinx.eu

## Instruction for Use

**PRODUCT NAME** Mercury-free glass thermometers

**INTENDED USE** The mercury-free glass thermometer has to be used for clinical medicine or body temperature measurement in the range of 35-42°C.

**APPLICABLE PARTS** oral, rectal, axillary

**CONTRAINDICATION** Can not be used for damaged skin.

### SPECIFICATION

Range: 35,0°C to 42,0°C.

Resolution: 0,1°C.

Accuracy: +0,1°C / -0,15°C.

Working Conditions: Temperature: 0°C to 42°C.

(It is only possible to measure temperatures that are above the ambient temperature!)

Storage Conditions: Temperature: 0°C to 42°C, Relative humidity ≤80 % RH.

### QUALITY GUARANTEE

- EU Council Directive 93/42/EEC MDD
- EN 12470-1:2000+A1:2009: Clinical thermometers - Part 1: Metallic liquid-in-glass thermometers with maximum device

### ⚠ PRECAUTIONS

1. Avoid falling and collision of the glass thermometer
2. Do not bend or bite the tip of the glass thermometer.
3. The glass thermometer shall be placed far away from children.
4. Children and disabled persons should be attended by an adult or medical staff during the entire temperature recording process.
5. The glass thermometer with a damaged outer casing must not be used, risk of injury.
6. The glass thermometer is only used to monitor body temperature, if find any abnormal body temperature, please visit your doctor immediately.

### INSTRUCTION FOR USE

1. Before measure, check the liquid column of the glass thermometer is below 36°C.
2. Clean the glass thermometer with 75% concentration alcohol before and after use.
3. Put the test port of the glass thermometer in the right part of the body (oral, axillary or rectal).
4. It needs 6 minutes for an accurate temperature measure, and then rotate the glass thermometer slowly back and forth to make the accurate reading. Within the measuring range, the measuring liquid column in the capillary tube shows the anthropometric temperature.
5. When the measure is completed, the measuring liquid must be returned to the bottom of the scale. In order to meet this requirement, it need to take the upper end of the thermometer as much as possible and throw the liquid column at least 5-12 times so as to reach below 36°C.

### DISPOSAL

If the glass thermometer is damaged and leaked, the spilled liquid can be removed by the brush, paper towel or gauze dipped in soapy water, which can be treated with the broken glass fragments for household waste.

### RELIABLE METHODS OF MEASURE TEMPERATURE

**Oral use:** measurement time: 6 minutes, normal temperature: approx. 37°C.

Doctors prefer oral measurement. It provides fast and accurate results. Place the thermometer probe to the left or right side under the tongue.

**Rectal use:** measurement time: 6 minutes, normal temperature: approx. 37,6°C.

Rectum measurement is preferred in the case of children. Insert the thermometer probe into the anus (approx. 2cm). You may use a little skin cream or baby oil on the head of the probe.

⚠ **STOP if you meet any resistance.**

**Axillary use:** measurement time: 6 minutes, normal temperature: approx. 36,7°C.

The axillary measurement method provides less accurate measurements than oral and rectal measurement. Wipe armpit with a dry towel. Place probe in the armpit and keep the arm pressed firmly against their side.

### NOTES ON CLEANING

1. The glass thermometer has a hermetically sealed, hygienic glass casing and can be disinfected without restriction.
2. The temperature of the cleaning and disinfectant solution must not exceed 42 °C.
3. Cold soapy water and 75% concentration alcohol can be used to clean the glass thermometer before and after use.
4. After cleaning, shall put it back to the rigid plastic pipes timely, so that it can be protected.
5. If this has already been used for rectal measurement, please mark this thermometer and separate storage. Shall not be used for oral use.

### STORAGE AND TRANSPORT

The device shall be stored at temperature 0 °C to 42 °C, relative humidity of less than 80%, the storage environment shall be kept non-corrosive gases.

### WARRANTY

Subject to normal use, the manufacturer grants an unlimited guarantee for the functionality and accuracy of this glass thermometer (except in the case of breakage of glass). Damage to the thermometer caused by improper treatment is not covered by the warranty.

Re-calibration is not necessary!

### SYMBOLS INSTRUCTION

	<b>CE mark</b>		<b>Lot number</b>
	<b>Manufacturer</b>		<b>European union representative</b>
	<b>Date of manufacture</b>		<b>Caution</b>

### CONTACT INFORMATION

 Jiaying Funing Medical Instrument Co.,Ltd.  
3<sup>rd</sup> Floor, Building C, Innovation town, Guigu Park,  
Luoxing Street, 314100 Jiashan County, Zhejiang Province, P.R. China

 Prolix GmbH  
Brehmstr. 56, 40239, Duesseldorf  
Tel.: 0049 2131 4051968-0  
Faks: 0049 2131 4051968-9  
Contact Person: G. Nianzhuang Liu  
Dimdi kod: DE/0000045300  
E-mail: dusai@eulinx.eu